

ADAM7 Antibody



PACO40558

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May play an important role in male reproduction including sperm maturation and gonadotrope function. This is a non catalytic metalloprotease-like protein.

Gene ID:

ADAM7

Uniprot

Q9H2U9

Synonyms:

Disintegrin and metalloproteinase domain-containing protein 7 (ADAM 7) (Sperm maturation-related glycoprotein GP-83), ADAM7, GP83


Immunogen:

Recombinant Human Disintegrin and metalloproteinase domain-containing protein 7 protein (177-210AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO40558 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human melanoma using PACO40558 at dilution of 1:100.