

# CCDC134 Antibody



PACO40578

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

**Gene ID:**

CCDC134

**Uniprot**

Q9H6E4

**Synonyms:**

Coiled-coil domain-containing protein 134, CCDC134

**Immunogen:**


Recombinant Human Coiled-coil domain-containing protein 134 protein (23-229AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO40578 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human brain tissue using PACO40578 at dilution of 1:100.