

# ZFP64 Antibody, FITC conjugated



PACO40636

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May be involved in transcriptional regulation.

**Gene ID:**

ZFP64

**Uniprot**

Q9NPA5

**Synonyms:**

Zinc finger protein 64 homolog, isoforms 1 and 2 (Zfp-64) (Zinc finger protein 338), ZFP64, ZNF338

**Immunogen:**

Recombinant Human Zinc finger protein 64 homolog, isoforms 1 and 2 protein (70-342AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A