

GJB4 Antibody, FITC conjugated



PACO40652

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

One gap junction consists of a cluster of closely packed pairs of transmembrane channels, the connexons, through which materials of low MW diffuse from one cell to a neighboring cell.

Gene ID:

GJB4

Uniprot

Q9NTQ9

Synonyms:

Gap junction beta-4 protein (Connexin-30.3) (Cx30.3), GJB4

Immunogen:

Recombinant Human Gap junction beta-4 protein (211-266AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A