

RBM12 Antibody, Biotin conjugated



PACO40657

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

nucleoplasm, poly(A) RNA binding.

Gene ID:

RBM12

Uniprot

Q9NTZ6

Synonyms:

RNA-binding protein 12 (RNA-binding motif protein 12) (SH3/WW domain anchor protein in the nucleus) (SWAN), RBM12, KIAA0765

Immunogen:

Recombinant Human RNA-binding protein 12 protein (1-200AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A