

C7orf50 Antibody, HRP conjugated



PACO40719

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	C7orf50
Source:	Uniprot
Rabbit	Q9BRJ6
Isotype:	Synonyms:
IgG	Uncharacterized protein C7orf50, C7orf50
Applications:	Immunogen:
ELISA	Recombinant Human Uncharacterized protein C7orf50 protein (1-194AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A