

# ZPBP Antibody, FITC conjugated



PACO40724

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May be implicated in gamete interaction during fertilization.

**Gene ID:**

ZPBP

**Uniprot**

Q9BS86

**Synonyms:**

Zona pellucida-binding protein 1 (Inner acrosomal membrane IAM38) (Sp38), ZPBP, ZPBP1

**Immunogen:**

Recombinant Human Zona pellucida-binding protein 1 protein (45-351AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A