

ZPBP Antibody, Biotin conjugated



PACO40725

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be implicated in gamete interaction during fertilization.

Gene ID:

ZPBP

Uniprot

Q9BS86

Synonyms:

Zona pellucida-binding protein 1 (Inner acrosomal membrane IAM38) (Sp38), ZPBP, ZPBP1

Immunogen:

Recombinant Human Zona pellucida-binding protein 1 protein (45-351AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A