

NICN1 Antibody, Biotin conjugated



PACO40729

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

NICN1

Uniprot

Q9BSH3

Synonyms:

Nicolin-1 (NPCEDRG) (Tubulin polyglutamylase complex subunit 5) (PGs5), NICN1

Immunogen:

Recombinant Human Nicolin-1 protein (1-213AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A