

PRRG3 Antibody



PACO40750

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

integral component of membrane.

Gene ID:

PRRG3

Uniprot

Q9BZD7

Synonyms:

Transmembrane gamma-carboxyglutamic acid, protein 3 (Proline-rich gamma-carboxyglutamic acid, protein 3) (Proline-rich Gla protein 3), PRRG3, PRGP3 TMG3


Immunogen:

Recombinant Human Transmembrane gamma-carboxyglutamic acid, protein 3 protein (102-231AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HeLa cells using PACO40750 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).