

# BRINP2 Antibody



PACO40758

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Inhibits neuronal cell proliferation by negative regulation of the cell cycle transition.

**Gene ID:**

BRINP2

**Uniprot**

Q9C0B6

**Synonyms:**

BMP/retinoic acid, inducible neural-specific protein 2 (DBCCR1-like protein 2), BRINP2, DBCCR1L2 FAM5B KIAA1747

**Immunogen:**


Recombinant Human BMP/retinoic acid, inducible neural-specific protein 2 protein (34-300AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO40758 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO40758 at dilution of 1:100.