

# ARV1 Antibody



PACO40790

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IF

**Recommended dilutions:**

ELISA:1:2000-1:10000, IF:1:50-1:200

**Protein Background:**

May act as a mediator of sterol homeostasis.

**Gene ID:**

ARV1

**Uniprot**

Q9H2C2

**Synonyms:**

Protein ARV1 (hARV1), ARV1

**Immunogen:**


Recombinant Human Protein ARV1 protein (1-143AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunofluorescent analysis of A2780 cells using PACO40790 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).