

TTYH1 Antibody



PACO40794

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IF:1:50-1:200

Protein Background:

Probable chloride channel. May be involved in cell adhesion.

Gene ID:

TTYH1

Uniprot

Q9H313

Synonyms:

Protein tweety homolog 1 (hTTY1), TTYH1


Immunogen:

Recombinant Human Protein tweety homolog 1 protein (50-178AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: TTYH1 antibody at 4 μ g/ml + Raji whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 50, 27, 36 kDa. Observed band size: 50 kDa.



Immunofluorescent analysis of HeLa cells using PACO40794 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).