

TRIT1 Antibody, HRP conjugated



PACO40799

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Catalyzes the transfer of a dimethylallyl group onto the adenine at position 37 of both cytosolic and mitochondrial tRNAs, leading to the formation of N6-(dimethylallyl)adenosine (i(6)A).

Gene ID:

TRIT1

Uniprot

Q9H3H1

Synonyms:

tRNA dimethylallyltransferase (EC 2.5.1.75) (Isopentenyl-diphosphate: tRNA isopentenyltransferase) (IPP transferase) (IPPT) (hGRO1) (tRNA isopentenyltransferase) (IPTase), TRIT1, IPT MOD5

Immunogen:

Recombinant Human tRNA dimethylallyltransferase protein (1-324AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A