

# TRIT1 Antibody, FITC conjugated



PACO40800

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Catalyzes the transfer of a dimethylallyl group onto the adenine at position 37 of both cytosolic and mitochondrial tRNAs, leading to the formation of N6-(dimethylallyl)adenosine (i(6)A).

**Gene ID:**

TRIT1

**Uniprot**

Q9H3H1

**Synonyms:**

tRNA dimethylallyltransferase (EC 2.5.1.75) (Isopentenyl-diphosphate: tRNA isopentenyltransferase) (IPP transferase) (IPPT) (hGRO1) (tRNA isopentenyltransferase) (IPTase), TRIT1, IPT MOD5

**Immunogen:**

Recombinant Human tRNA dimethylallyltransferase protein (1-324AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A