

SLITRK6 Antibody, Biotin conjugated



PACO40809

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Regulator of neurite outgrowth required for normal hearing and vision.

Gene ID:

SLITRK6

Uniprot

Q9H5Y7

Synonyms:

SLIT and NTRK-like protein 6, SLITRK6

Immunogen:

Recombinant Human SLIT and NTRK-like protein 6 protein (621-821AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A