

NSUN3 Antibody, FITC conjugated



PACO40812

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May have S-adenosyl-L-methionine-dependent methyl-transferase activity.

Gene ID:

NSUN3

Uniprot

Q9H649

Synonyms:

tRNA (cytosine(34)-C(5))-methyltransferase, mitochondrial (EC 2.1.1. -)
(NOL1/NOP2/Sun domain family member 3), NSUN3

Immunogen:

Recombinant Human tRNA (cytosine(34)-C(5))-methyltransferase, mitochondrial protein (1-200AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A