

C12orf49 Antibody



PACO40814

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

C12orf49

Uniprot

Q9H741

Synonyms:

UPF0454 protein C12orf49, C12orf49


Immunogen:

Recombinant Human UPF0454 protein C12orf49 protein (34-205AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO40814 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO40814 at dilution of 1:100.