

RHOF Antibody, Biotin conjugated



PACO40833

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plasma membrane-associated small GTPase which cycles between an active GTP-bound and an inactive GDP-bound state. Causes the formation of thin, actin-rich surface projections called filopodia. Functions cooperatively with CDC42 and Rac to generate additional structures, increasing the diversity of actin-based morphology.

Gene ID:

RHOF

Uniprot

Q9HBH0

Synonyms:

Rho-related GTP-binding protein RhoF (Rho family GTPase Rif) (Rho in filopodia), RHOF, ARHF RIF

Immunogen:

Recombinant Human Rho-related GTP-binding protein RhoF protein (1-170AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A