

PLAC1 Antibody



PACO40834

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in placental development.

Gene ID:

PLAC1

Uniprot

Q9HBJ0

Synonyms:

Placenta-specific protein 1, PLAC1

Immunogen:

Recombinant Human Placenta-specific protein 1 protein (23-212AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A