

NRN1 Antibody, HRP conjugated



PACO40847

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Promotes neurite outgrowth and especially branching of neuritic processes in primary hippocampal and cortical cells.

Gene ID:

NRN1

Uniprot

Q9NPD7

Synonyms:

Neuritin, NRN1, NRN

Immunogen:

Recombinant Human Neuritin protein (28-116AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A