

# NRN1 Antibody, Biotin conjugated



PACO40849

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Promotes neurite outgrowth and especially branching of neuritic processes in primary hippocampal and cortical cells.

**Gene ID:**

NRN1

**Uniprot**

Q9NPD7

**Synonyms:**

Neuritin, NRN1, NRN

**Immunogen:**

Recombinant Human Neuritin protein (28-116AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A