

NRIP3 Antibody



PACO40854

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

NRIP3

Uniprot

Q9NQ35

Synonyms:

Nuclear receptor-interacting protein 3 (Sarcoma antigen NY-SAR-105), NRIP3, C11orf14


Immunogen:

Recombinant Human Nuclear receptor-interacting protein 3 protein (1-241AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human liver cancer using PACO40854 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO40854 at dilution of 1:100.