

PPAN Antibody, FITC conjugated



PACO40860

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May have a role in cell growth.

Gene ID:

PPAN

Uniprot

Q9NQ55

Synonyms:

Suppressor of SWI4 1 homolog (Ssf-1) (Brix domain-containing protein 3) (Peter Pan homolog), PPAN, BXDC3 SSF1

Immunogen:

Recombinant Human Suppressor of SWI4 1 homolog protein (1-473AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A