

OXSM Antibody, Biotin conjugated



PACO40893

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May play a role in the biosynthesis of lipoic acid, as well as longer chain fatty acid, required for optimal mitochondrial function.

Gene ID:

OXSM

Uniprot

Q9NWU1

Synonyms:

3-oxoacyl-[acyl-carrier-protein] synthase, mitochondrial (EC 2.3.1.41) (Beta-ketoacyl-ACP synthase), OXSM

Immunogen:

Recombinant Human 3-oxoacyl-[acyl-carrier-protein] synthase, mitochondrial protein (78-345AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A