

GDAP2 Antibody



PACO40902

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:

Gene ID:

GDAP2

Uniprot

Q9NXN4

Synonyms:

Ganglioside-induced differentiation-associated protein 2, GDAP2


Immunogen:

Recombinant Human Ganglioside-induced differentiation-associated protein 2 protein (1-496AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: GDAP2 antibody at 2 μ g/ml + MCF-7 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 57 kDa. Observed band size: 57 kDa.