

GDAP2 Antibody, Biotin conjugated



PACO40905

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

GDAP2

Uniprot

Q9NXN4

Synonyms:

Ganglioside-induced differentiation-associated protein 2, GDAP2

Immunogen:

Recombinant Human Ganglioside-induced differentiation-associated protein 2 protein (1-496AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A