

UTP6 Antibody, HRP conjugated



PACO40911

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in nucleolar processing of pre-18S ribosomal RNA.

Gene ID:

UTP6

Uniprot

Q9NYH9

Synonyms:

U3 small nucleolar RNA-associated protein 6 homolog (Hepatocellular carcinoma-associated antigen 66) (Multiple hat domains protein), UTP6, C17orf40 HCA66 MHAT

Immunogen:

Recombinant Human U3 small nucleolar RNA-associated protein 6 homolog protein (1-220AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A