

TAS2R14 Antibody, HRP conjugated



PACO40915

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Receptor that may play a role in the perception of bitterness and is gustducin-linked. May play a role in sensing the chemical composition of the gastrointestinal content. The activity of this receptor may stimulate alpha gustducin, mediate PLC-beta-2 activation and lead to the gating of TRPM5.

Gene ID:

TAS2R14

Uniprot

Q9NYV8

Synonyms:

Taste receptor type 2 member 14 (T2R14) (Taste receptor family B member 1) (TRB1), TAS2R14

Immunogen:

Recombinant Human Taste receptor type 2 member 14 protein (283-317AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A