

MBD5 Antibody, FITC conjugated



PACO40932

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Binds to heterochromatin. Does not interact with either methylated or unmethylated DNA (in vitro).

Gene ID:

MBD5

Uniprot

Q9P267

Synonyms:

Methyl-CpG-binding domain protein 5 (Methyl-CpG-binding protein MBD5), MBD5, KIAA1461

Immunogen:

Recombinant Human Methyl-CpG-binding domain protein 5 protein (623-840AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A