

# CEP126 Antibody, Biotin conjugated



PACO40937

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May participate in cytokinesis.

**Gene ID:**

CEP126

**Uniprot**

Q9P2H0

**Synonyms:**

Centrosomal protein of 126 kDa, CEP126, KIAA1377

**Immunogen:**

Recombinant Human Centrosomal protein of 126 kDa protein (559-670AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A