

DCTN5 Antibody, FITC conjugated



PACO40964

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

centrosome, cytosol, nuclear membrane, nucleoplasm, antigen processing and presentation of exogenous peptide antigen via MHC class II.

Gene ID:

DCTN5

Uniprot

Q9BTE1

Synonyms:

Dynactin subunit 5 (Dynactin subunit p25), DCTN5

Immunogen:

Recombinant Human Dynactin subunit 5 protein (1-182AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A