

RTL8A Antibody, FITC conjugated



PACO40984

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

RTL8A

Uniprot

Q9BWD3

Synonyms:

Retrotransposon Gag-like protein 8A (Mammalian retrotransposon derived protein 8A),
RTL8A, CXX1B FAM127B MAR8A

Immunogen:

Recombinant Human Retrotransposon Gag-like protein 8A protein (1-113AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A