

# RHPN1-AS1 Antibody



PACO40986

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

**Gene ID:**

RHPN1-AS1

**Uniprot**

Q9BWJ2

**Synonyms:**

Putative uncharacterized protein encoded by RHPN1-AS1 (RHPN1 antisense RNA 1) (RHPN1 antisense gene protein 1), RHPN1-AS1, C8orf51

**Immunogen:**


Recombinant Human Putative uncharacterized protein encoded by RHPN1-AS1 protein (1-59AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human brain tissue using PACO40986 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO40986 at dilution of 1:100.