

RHPN1-AS1 Antibody, FITC conjugated



PACO40988

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

RHPN1-AS1

Uniprot

Q9BWJ2

Synonyms:

Putative uncharacterized protein encoded by RHPN1-AS1 (RHPN1 antisense RNA 1) (RHPN1 antisense gene protein 1), RHPN1-AS1, C8orf51

Immunogen:

Recombinant Human Putative uncharacterized protein encoded by RHPN1-AS1 protein (1-59AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A