

# C1QTNF1 Antibody, FITC conjugated



PACO40996

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

C1QTNF1

**Uniprot**

Q9BXJ1

**Synonyms:**

Complement C1q tumor necrosis factor-related protein 1 (G protein-coupled receptor-interacting protein) (GIP), C1QTNF1, CTRP1

**Immunogen:**

Recombinant Human Complement C1q tumor necrosis factor-related protein 1 protein (26-281AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A