

# NLRP1 Antibody, FITC conjugated



PACO41028

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NLRP1

**Uniprot**

Q9C000

**Synonyms:**

NACHT, LRR and PYD domains-containing protein 1 (Caspase recruitment domain-containing protein 7) (Death effector filament-forming ced-4-like apoptosis protein) (Nucleotide-binding domain and caspase recruitment domain), NLRP1, CARD7 DEFCAP KIAA0926 NAC NALP1

**Immunogen:**

Recombinant Human NACHT, LRR and PYD domains-containing protein 1 protein (1-146AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A