

SETD5 Antibody, FITC conjugated



PACO41032

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

SETD5

Uniprot

Q9C0A6

Synonyms:

SET domain-containing protein 5, SETD5, KIAA1757

Immunogen:

Recombinant Human SET domain-containing protein 5 protein (519-700AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A