

SDCBP2 Antibody, HRP conjugated



PACO41055

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cytoplasm, extracellular exosome, plasma membrane, protein C-terminus binding, protein heterodimerization activity, protein homodimerization activity, intracellular signal transduction, intracellular transport, nervous system development.

Gene ID:

SDCBP2

Uniprot

Q9H190

Synonyms:

Syntenin-2 (Similar to TACIP18) (SITAC) (Syndecan-binding protein 2), SDCBP2, SITAC18

Immunogen:

Recombinant Human Syntenin-2 protein (11-207AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A