

# R3HDML Antibody



PACO41066

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Putative serine protease inhibitor.

**Gene ID:**

R3HDML

**Uniprot**

Q9H3Y0

**Synonyms:**

Peptidase inhibitor R3HDML (Cysteine-rich secretory protein R3HDML), R3HDML

**Immunogen:**

Recombinant Human Peptidase inhibitor R3HDML protein (90-253AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A