

ACTR8 Antibody, HRP conjugated



PACO41111

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays an important role in the functional organization of mitotic chromosomes. Exhibits low basal ATPase activity, and unable to polymerize.

Gene ID:

ACTR8

Uniprot

Q9H981

Synonyms:

Actin-related protein 8 (hArap8) (INO80 complex subunit N), ACTR8, ARP8 INO80N

Immunogen:

Recombinant Human Actin-related protein 8 protein (35-175AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A