

ANKRA2 Antibody



PACO41114

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May facilitate endocytosis by linking megalin to components of the cytoskeleton or endocytic machinery.

Gene ID:

ANKRA2

Uniprot

Q9H9E1

Synonyms:

Ankyrin repeat family A protein 2 (RFXANK-like protein 2), ANKRA2, ANKRA


Immunogen:

Recombinant Human Ankyrin repeat family A protein 2 protein (1-313AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human liver tissue using PACO41114 at dilution of 1:100.