

ZDHHC4 Antibody, FITC conjugated



PACO41156

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

endoplasmic reticulum, Golgi apparatus, protein-cysteine S-palmitoyltransferase activity, peptidyl-L-cysteine S-palmitoylation, protein targeting to membrane.

Gene ID:

ZDHHC4

Uniprot

Q9NPG8

Synonyms:

Probable palmitoyltransferase ZDHHC4 (EC 2.3.1.225) (Zinc finger DHHC domain-containing protein 4) (DHHC-4) (Zinc finger protein 374), ZDHHC4, ZNF374

Immunogen:

Recombinant Human Probable palmitoyltransferase ZDHHC4 protein (277-344AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A