

# CYTL1 Antibody, FITC conjugated



PACO41168

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

extracellular space, receptor binding, signal transduction.

**Gene ID:**

CYTL1

**Uniprot**

Q9NRR1

**Synonyms:**

Cytokine-like protein 1 (Protein C17), CYTL1, C4orf4

**Immunogen:**

Recombinant Human Cytokine-like protein 1 protein (24-136AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A