

# OBP2A Antibody



PACO41202

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Probably binds and transports small hydrophobic volatile molecules with a higher affinity for aldehydes and large fatty acid, .

**Gene ID:**

OBP2A

**Uniprot**

Q9NY56

**Synonyms:**

Odorant-binding protein 2a (Odorant-binding protein IIa) (OBPIIa), OBP2A

**Immunogen:**


Recombinant Human Odorant-binding protein 2a protein (16-170AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human liver cancer using PACO41202 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human pancreatic cancer using PACO41202 at dilution of 1:100.