

COL20A1 Antibody, FITC conjugated



PACO41220

Product Information

Size:	Protein Background:
50ug	Probable collagen protein.
Reactivity:	Gene ID:
Human	COL20A1
Source:	Uniprot
Rabbit	Q9P218
Isotype:	Synonyms:
IgG	Collagen alpha-1(XX) chain, COL20A1, KIAA1510
Applications:	Immunogen:
ELISA	Recombinant Human Collagen alpha-1(XX) chain protein (1143-1297AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A