

WRKY40 Antibody



PACO41230

Product Information

Size:

50ug

Reactivity:

Arabidopsis thaliana

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Transcription factor. Interacts specifically with the W box (5'-(T)TGAC[CT]-3'), a frequently occurring elicitor-responsive cis-acting element.

Gene ID:

WRKY40

Uniprot

Q9SAH7

Synonyms:

Probable WRKY transcription factor 40 (WRKY DNA-binding protein 40), WRKY40

Immunogen:

Recombinant Arabidopsis thaliana Probable WRKY transcription factor 40 protein (1-302AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A