

ABHD14A Antibody, Biotin conjugated



PACO41253

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Possible role in granule neuron development.

Gene ID:

ABHD14A

Uniprot

Q9BUJ0

Synonyms:

Protein ABHD14A (EC 3) (Alpha/beta hydrolase domain-containing protein 14A) (Abhydrolase domain-containing protein 14A), ABHD14A

Immunogen:

Recombinant Human Alpha/β hydrolase domain-containing protein 14A protein (56-258AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A