

MRPL1 Antibody, HRP conjugated



PACO41275

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

MRPL1

Uniprot

Q9BYD6

Synonyms:

39S ribosomal protein L1, mitochondrial (L1mt) (MRP-L1) (Mitochondrial large ribosomal subunit protein uL1m), MRPL1

Immunogen:

Recombinant Human 39S ribosomal protein L1, mitochondrial protein (51-324AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A