

# MRPL1 Antibody, FITC conjugated



PACO41276

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MRPL1

**Uniprot**

Q9BYD6

**Synonyms:**

39S ribosomal protein L1, mitochondrial (L1mt) (MRP-L1) (Mitochondrial large ribosomal subunit protein uL1m), MRPL1

**Immunogen:**

Recombinant Human 39S ribosomal protein L1, mitochondrial protein (51-324AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A